

# HQ 90



## TOUCHLESS WATER PURIFICATION SYSTEM



### TOUCHLESS SENSOR ACTIVATION

Hygienic, touch-free operation.



### ACTIVATED OXYGEN SELF-CLEANING

Advanced technology circulates water and cleans tanks and lines to maintain water freshness.



### ANTI-MICROBIAL SURFACES

Anti-microbial additive maintains cleanliness of unit exterior.




### 24/7 SYSTEM MONITORING

Continuous monitoring of system performance and water quality with real-time status on LCD screen.



### ADVANCED PURIFICATION TECHNOLOGY

Proprietary, multi-stage, reverse osmosis process is certified to industry standards.



### HIGH-CAPACITY DRIP TRAY

Optional drain connection prevents need for emptying.



### ENERGY-SAVING MODE

Eco-Mode for reduced power consumption.



Hydrate HQ of Mid Atlantic  
1550 Caton Center Dr. Ste. D  
Baltimore, MD 21227  
443.873.0300

Hydrate HQ of SoCal & San Diego  
9340 Santa Anita Ave. Ste. 103  
Rancho Cucamonga, CA 91730  
909.454.3400

# HQ90

## TOUCHLESS WATER PURIFICATION SYSTEM

Large LCD screen displays real-time system performance and water quality status.

High-end glass front panel

Touch-free sensor activation


Hot water safety lock

Color-changing indicator light for precise dispensing

Large drip tray

High-capacity cold tank

Dent-resistant steel frame



## Specifications

Weight	69.9 lbs (31.7 kgs)
Dimensions	14" W x 50.9" H x 17" L (355.6mm W x 1293mm H x 431.8mm L)
Cold Tank Capacity	3 gallons (11 liters)
Hot Tank Capacity	0.5 gallons (2 liters)
Dispense Area	11.5" tall (293mm)
Recommended Water Pressure	40-80 PSI
Rated Voltage/Frequency	120v/60hz   220v/50hz
Power Consumption	Hot: 2.5A/300W   Cold: 1.0A/120W

## The Hydrate HQ Purification Process



Hydrate HQ of Mid Atlantic  
1550 Caton Center Dr. Ste. D  
Baltimore, MD 21227  
443.873.0300

Hydrate HQ of SoCal & San Diego  
9340 Santa Anita Ave. Ste. 103  
Rancho Cucamonga, CA 91730  
909.454.3400